

Pulse Programming Bootcamp

Day 2: Pulse Sequences

Mike Tyszka, Ph.D.



Pulse Programming Bootcamp 2005


Caltech Brain Imaging Center

California Institute of Technology

2D Multislice vs 3D


Gradient Echo Sequences

Full Gradient Echo Pulse Sequence



Steady State Free Precession

Longitudinal Steady-State



RF pulse rotates M_z magnetization about x-axis

$$M_z(0^+) = M_z(0^-) \cos(\alpha)$$


T_1 relaxation during free-precession period (TR)

$$M_z(TR) = M_{xy}(0^+) e^{-TR/T_1} + M_0(1 - e^{-TR/T_1})$$

Equate $M_z(0^-)$ with $M_z(TR)$ in steady-state assuming complete M_{xy} spoiling

$$M_z^+ = M_0 \frac{1 - e^{-TR/T_1}}{1 - \cos(\alpha) e^{-TR/T_1}}$$

The Ernst Angle



$$\theta_E = \cos^{-1} \left(e^{-TR/T_1} \right)$$


Maximum signal
obtained at Ernst
angle, **not** 90°
when $TR < \text{or } \approx T_I$

$$M_{xy}^+ = M_0 \frac{1 - e^{-TR/T_1}}{1 - \cos(\alpha) e^{-TR/T_1}} \sin(\alpha)$$


Full Steady-State

Transverse Magnetization Preserved. $TR < T_2$


Both M_z and M_{xy} attain a steady-state.




Steady-State Gradient Echo Sequence




TrueFISP Pulse Sequence



Off-resonance Effects in Full Steady-State




TrueFISP Isofield Null Points




Spin Echo Sequences

The Hahn Spin Echo Explained




Spin Echo Pulse Sequence




RARE and Turbo Spin Echo

Spin Echo Pulse Sequence




Carr-Purcell-Meiboom-Gill (CPMG) Sequence




Multiple spin echoes refocused
from single excitation


Turbo Spin Echo Pulse Sequence



TSE k-space trajectory




TSE k-space trajectory (continued)




Stimulated Echoes

Stimulated Echo Sequence





Stimulated Echo Imaging




Echo Planar Imaging

Echo-planar Imaging



Gradient Echo EPI




GE-EPI Images


Gradient Echo EPI Images

GE-EPI



Shots : 1
TR : 3000ms
TE : 60ms
Slice : 5mm/2.5mm (18)
Matrix : 128 x 128
FOV : 24cm x 24cm
Time : 6s
NEX : 1



Spin Echo EPI




SE-EPI Images



SE-EPI


Shots : 1
TR : 3000ms
TE : 60ms
Slice : 5mm/2.5mm (18)
Matrix : 128 x 128
FOV : 24cm x 24cm
Time : 3s
NEX : 1

EPI Bandwidths



Static Inhomogeneity and EPI

Frequency Encoding Direction



Phase Encoding Direction


Voxel displacement
due to static field
inhomogeneity

$$\Delta y = \frac{\gamma \Delta B_0 \cdot \text{FOV}_y}{\text{BW}_y}$$


Off-resonance Frequency (Hz) at 1.5T




Effect of Echo Train Length



Ramp Sampling




Without Ramp Sampling




With Ramp Sampling

Multi-shot EPI k-space Trajectory


Single Shot



Two Shot



Effect of Multi-shot




Increasing the number of shots per image decreases the EPI echo-train length per shot.

SE-EPI


Shots	: 1-8
TR	: 3000ms
TE	: 60ms
Slice	: 5mm/2.5mm (18)
Matrix	: 256 x 256
FOV	: 24cm x 24cm
Time	: 12s-27s
NEX	: 1

Magnetization Preparation




Longitudinal or transverse magnetization carries magnetization contrast prepared by initial pulse sequence

Inversion-Recovery EPI




IR-EPI Images




IR-SE-EPI

Shots : 1
TI : 400ms
TR : 3000ms
TE : 60ms
Slice : 5mm/2.5mm (18)
Matrix : 128 x 128
FOV : 24cm x 24cm
Time : 12s
NEX : 1

MP-RAGE Pulse Sequence




T₁-weighted MP-RAGE




Phase Encode Ordering


Sequential



Interleaved





Centric





Segmented Acquisitions


Segmented Acquisition (eg TSE/RARE)



Echo 4 

Echo 3 

Echo 2 

Echo 1 

Segment Ordering?

Intrasegment Ordering?

Sequential

Sequential

**Tomorrow :
Corrections and Compensations**